

Anti-CD235 antibody [HIR2] {PE/Cyanine5} (STJA0001522)

STJA0001522

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD235 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	HIR2
Concentration	5 µg L/Test
Conjugation	PE/Cyanine5
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2bk
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	2994
	2993
Gene Symbol	GYPB
	GYPA
Uniprot ID	GLPB_HUMAN
	GLPA_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	