

## Anti-CD11b antibody [M1/70] {PE/Cyanine5.5} (STJA0001526)

STJA0001526

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rat monoclonal antibody anti-CD11b is suitable for use in Flow Cytometry research applications.
<b>Applications</b>	FC
<b>Host/Source</b>	Rat
<b>Reactivity</b>	Human/Mouse

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	M1/70
<b>Concentration</b>	5 µg L/Test
<b>Conjugation</b>	PE/Cyanine5.5
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
<b>Isotype</b>	IgG2bk
<b>Storage Instruction</b>	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

### TARGET INFORMATION

<b>Gene ID</b>	NA <a href="#">3684</a>
<b>Gene Symbol</b>	<a href="#">Itgam</a> <a href="#">ITGAM</a>
<b>Uniprot ID</b>	<a href="#">ITAM_MOUSE</a> <a href="#">ITAM_HUMAN</a>
<b>Immunogen</b>	
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	