

Anti-CD64 antibody [10.1] (STJA0001528)

STJA0001528

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD64 is suitable for use in Flow Cytometry and Immunohistochemistry research applications.
Applications	FC/IHC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	10.1
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300
Isotype	IgG1k
Storage Instruction	Store at 2-8°C. Do not freeze this product.

TARGET INFORMATION

Gene ID	NA 2209
Gene Symbol	FCGR1BP FCGR1A
Uniprot ID	FCGRB_HUMAN FCGR1_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	