

Anti-CD103/Itgae antibody [M290] {Biotin} (STJA0001552)

STJA0001552

GENERAL INFORMATION

| | |
|--------------------------|------------------------------------------------------------------------------------------------------------|
| Product Type | Primary antibodies |
| Short Description | Rat monoclonal antibody anti-Integrin alpha-E is suitable for use in Flow Cytometry research applications. |
| Applications | FC |
| Host/Source | Rat |
| Reactivity | Mouse |

PRODUCT PROPERTIES

| | |
|----------------------------|---------------------------------------------------------------------------------------------|
| Clonality | Monoclonal |
| Clone ID | M290 |
| Concentration | 0.5 mg/mL |
| Conjugation | Biotin |
| Purification | |
| Dilution Range | |
| Formulation | Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant. |
| Isotype | IgG2ak |
| Storage Instruction | Store at 2-8°C. Do not freeze this product. |

TARGET INFORMATION

| | |
|---------------------------|----------------------------|
| Gene ID | 16407 |
| Gene Symbol | Itgae |
| Uniprot ID | ITAE_MOUSE |
| Immunogen | |
| Immunogen Region | |
| Specificity | |
| Immunogen Sequence | |