

**Anti-CD106/Vcam1 antibody [M/K-2.7] (STJA0001564)**

STJA0001564

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rat monoclonal antibody anti-Vascular cell adhesion protein 1 is suitable for use in Flow Cytometry and Immunofluorescence research applications.
<b>Applications</b>	FC/IF
<b>Host/Source</b>	Rat
<b>Reactivity</b>	Mouse

**PRODUCT PROPERTIES**

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	M/K-2.7
<b>Concentration</b>	0.5 mg/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	PBS with 0.05% Proclin300
<b>Isotype</b>	IgG1k
<b>Storage</b>	Store at 2-8°C. Do not freeze this product.
<b>Instruction</b>	

**TARGET INFORMATION**

<b>Gene ID</b>	
<b>Gene Symbol</b>	<a href="#">Vcam1</a>
<b>Uniprot ID</b>	<a href="#">VCAM1_MOUSE</a>
<b>Immunogen</b>	
<b>Immunogen</b>	
<b>Region</b>	
<b>Specificity</b>	
<b>Immunogen</b>	
<b>Sequence</b>	