

**Anti-CD106/Vcam1 antibody [M/K-2.7] {PE/Cyanine5.5} (STJA0001569)**

STJA0001569

**GENERAL INFORMATION**

|                          |                                                                                                                            |
|--------------------------|----------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type</b>      | Primary antibodies                                                                                                         |
| <b>Short Description</b> | Rat monoclonal antibody anti-Vascular cell adhesion protein 1 is suitable for use in Flow Cytometry research applications. |
| <b>Applications</b>      | FC                                                                                                                         |
| <b>Host/Source</b>       | Rat                                                                                                                        |
| <b>Reactivity</b>        | Mouse                                                                                                                      |

**PRODUCT PROPERTIES**

|                            |                                                                                                              |
|----------------------------|--------------------------------------------------------------------------------------------------------------|
| <b>Clonality</b>           | Monoclonal                                                                                                   |
| <b>Clone ID</b>            | M/K-2.7                                                                                                      |
| <b>Concentration</b>       | 0.2 mg/mL                                                                                                    |
| <b>Conjugation</b>         | PE/Cyanine5.5                                                                                                |
| <b>Purification</b>        |                                                                                                              |
| <b>Dilution Range</b>      |                                                                                                              |
| <b>Formulation</b>         | Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.                  |
| <b>Isotype</b>             | IgG1k                                                                                                        |
| <b>Storage Instruction</b> | Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product. |

**TARGET INFORMATION**

|                           |                             |
|---------------------------|-----------------------------|
| <b>Gene ID</b>            |                             |
| <b>Gene Symbol</b>        | <a href="#">Vcam1</a>       |
| <b>Uniprot ID</b>         | <a href="#">VCAM1_MOUSE</a> |
| <b>Immunogen</b>          |                             |
| <b>Immunogen Region</b>   |                             |
| <b>Specificity</b>        |                             |
| <b>Immunogen Sequence</b> |                             |