

Anti-Cd4 antibody [GK1.5] (STJA0001615)

STJA0001615

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-T-cell surface glycoprotein CD4 is suitable for use in Universal, Universal, Flow Cytometry, Western Blot, Immunohistochemistry and Immunoprecipitation research applications.
Applications	Depletion/Block/FC/WB/IHC/IP
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	GK1.5
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	0.2 µm filtered in PBS, pH7.2. Azide Free (AF) /Low Endotoxin (LE) : Contains no stabilizers or preservatives. Endotoxin level is < 2 EU/mg as determined by LAL gel clotting assay.
Isotype	IgG2bk
Storage	Store at 2-8°C. Do not freeze this product.
Instruction	

TARGET INFORMATION

Gene ID	12504
Gene Symbol	Cd4
Uniprot ID	CD4_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	