

Anti-CD44 antibody [IM7] (STJA0001633)

STJA0001633

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-CD44 is suitable for use in Universal, Flow Cytometry, Immunohistochemistry and Immunoprecipitation research applications.
Applications	Neut/FC/IHC/IP
Host/Source	Rat
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	IM7
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	0.2 µm filtered in PBS, pH7.2. Azide Free (AF) /Low Endotoxin (LE) : Contains no stabilizers or preservatives. Endotoxin level is < 2 EU/mg as determined by LAL gel clotting assay.
Isotype	IgG2bk
Storage	Store at 2-8°C. Do not freeze this product.
Instruction	

TARGET INFORMATION

Gene ID	960
	12505
Gene Symbol	CD44
	Cd44
Uniprot ID	CD44_HUMAN
	CD44_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	