

Anti-CD25/Il2ra antibody [PC-61.5.3] (STJA0001654)

STJA0001654

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Interleukin-2 receptor subunit alpha is suitable for use in Flow Cytometry, Immunohistochemistry and Immunoprecipitation research applications.
Applications	FC/IHC/IP
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	PC-61.5.3
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300
Isotype	IgG1k
Storage	Store at 2-8°C. Do not freeze this product.
Instruction	

TARGET INFORMATION

Gene ID	16184
Gene Symbol	Il2ra
Uniprot ID	IL2RA_MOUSE
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	