

**Anti-CD3 Epsilon antibody [145-2C11] {PE/Cyanine5.5} (STJA0001668)**

STJA0001668

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Armenian Hamster monoclonal antibody anti-T-cell surface glycoprotein CD3 epsilon chain is suitable for use in Flow Cytometry research applications.
<b>Applications</b>	FC
<b>Host/Source</b>	Armenian Hamster
<b>Reactivity</b>	Mouse

**PRODUCT PROPERTIES**

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	145-2C11
<b>Concentration</b>	0.2 mg/mL
<b>Conjugation</b>	PE/Cyanine5.5
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
<b>Isotype</b>	IgG
<b>Storage</b>	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.
<b>Instruction</b>	

**TARGET INFORMATION**

<b>Gene ID</b>	<a href="#">12501</a>
<b>Gene Symbol</b>	<a href="#">Cd3e</a>
<b>Uniprot ID</b>	<a href="#">CD3E_MOUSE</a>
<b>Immunogen</b>	
<b>Immunogen</b>	
<b>Region</b>	
<b>Specificity</b>	
<b>Immunogen</b>	
<b>Sequence</b>	