

Anti-Cd8a antibody [53-6.7] (STJA0001670)

STJA0001670

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-T-cell surface glycoprotein CD8 alpha chain is suitable for use in Universal, Universal, Flow Cytometry, Western Blot, Immunofluorescence and Immunoprecipitation research applications.
Applications	Block/Depletion/FC/WB/IF/IP
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	53-6.7
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	
Dilution	
Range	
Formulation	0.2 µm filtered in PBS, pH7.2. Azide Free (AF) /Low Endotoxin (LE) : Contains no stabilizers or preservatives. Endotoxin level is < 2 EU/mg as determined by LAL gel clotting assay.
Isotype	IgG2ak
Storage	Store at 2-8°C. Do not freeze this product.
Instruction	

TARGET INFORMATION

Gene ID	12525
Gene Symbol	Cd8a
Uniprot ID	CD8A_MOUSE
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	