

Anti-CD115/Csf1r antibody [AFS98] {APC} (STJA0001692)

STJA0001692

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Macrophage colony-stimulating factor 1 receptor is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	AFS98
Concentration	5 Mu L/Test
Conjugation	APC
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2ak
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	12978
Gene Symbol	Csf1r
Uniprot ID	CSF1R_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen	
Sequence	