

Anti-CD4 antibody [RPA-T4] (STJA0001700)

STJA0001700

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-T-cell surface glycoprotein CD4 is suitable for use in Universal, Flow Cytometry, Immunohistochemistry and Immunofluorescence research applications.
Applications	Block/FC/IHC/IF
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	RPA-T4
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	0.2 µm filtered in PBS, pH7.2. Azide Free (AF) /Low Endotoxin (LE) : Contains no stabilizers or preservatives. Endotoxin level is < 2 EU/mg as determined by LAL gel clotting assay.
Isotype	IgG1k
Storage	Store at 2-8°C. Do not freeze this product.
Instruction	

TARGET INFORMATION

Gene ID	920
Gene Symbol	CD4
Uniprot ID	CD4_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	