

Anti-CD45R/B220 antibody [RA3.3A 1/6.1] (STJA0001711)

STJA0001711

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-CD45R/B220 is suitable for use in Universal, Universal, Flow Cytometry, Immunohistochemistry and Immunoprecipitation research applications.
Applications	Cell Sep-Neg/Depletion/FC/IHC/IP
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	RA3.3A 1/6.1
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	0.2 µm filtered in PBS, pH7.2. Azide Free (AF) /Low Endotoxin (LE) : Contains no stabilizers or preservatives. Endotoxin level is < 2 EU/mg as determined by LAL gel clotting assay.
Isotype	IgMk
Storage	Store at 2-8°C. Do not freeze this product.
Instruction	

TARGET INFORMATION

Gene ID	19264 5788
Gene Symbol	Ptprc PTPRC
Uniprot ID	PTPRC_MOUSE PTPRC_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	