

Anti-CD45R/B220 antibody [RA3.3A 1/6.1] {Biotin} (STJA0001713)

STJA0001713

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-CD45R/B220 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	RA3.3A 1/6.1
Concentration	0.5 mg/mL
Conjugation	Biotin
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgMk
Storage Instruction	Store at 2-8°C. Do not freeze this product.

TARGET INFORMATION

Gene ID	19264 5788
Gene Symbol	Ptprc PTPRC
Uniprot ID	PTPRC_MOUSE PTPRC_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	