

Anti-CD18/Itgb2 antibody [M18/2] {PE} (STJA0001724)

STJA0001724

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Integrin beta-2 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	M18/2
Concentration	0.2 mg/mL
Conjugation	PE
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2ak
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	16414
Gene Symbol	Itgb2
Uniprot ID	ITB2_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	