

Anti-CD183/Cxcr3 antibody [CXCR3-173] {PE/Cyanine7} (STJA0001733)

STJA0001733

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Armenian Hamster monoclonal antibody anti-C-X-C chemokine receptor type 3 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Armenian Hamster
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	CXCR3-173
Concentration	0.2 mg/mL
Conjugation	PE/Cyanine7
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	12766
Gene Symbol	Cxcr3
Uniprot ID	CXCR3_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen	
Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081