

Anti-CD119/Ifngr1 antibody [GR-20] {PE} (STJA0001739)
STJA0001739**GENERAL INFORMATION**

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Interferon gamma receptor 1 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	GR-20
Concentration	5 Mu L/Test
Conjugation	PE
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2ak
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	15979
Gene Symbol	Ifngr1
Uniprot ID	INGR1_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	