

Anti-Ly-6G/Ly-6C (Gr-1) antibody [RB6-8C5] (STJA0001754)

STJA0001754

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Ly-6G/Ly-6C (Gr-1) is suitable for use in Flow Cytometry, Western Blot, Immunohistochemistry and Immunoprecipitation research applications.
Applications	FC/WB/IHC/IP
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	RB6-8C5
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300
Isotype	IgG2bk
Storage	Store at 2-8°C. Do not freeze this product.
Instruction	

TARGET INFORMATION

Gene ID	546644 17067 100041546
Gene Symbol	Ly6g Ly6c1 LY6G_MOUSE LY6C1_MOUSE LY6C2_MOUSE <
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	