

Anti-Ly6C antibody [Monts 1] (STJA0001758)

STJA0001758

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Ly6C is suitable for use in Universal and Flow Cytometry research applications.
Applications	Depletion/FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	Monts 1
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	0.2 µm filtered in PBS, pH7.2. Azide Free (AF) /Low Endotoxin (LE) : Contains no stabilizers or preservatives. Endotoxin level is < 2 EU/mg as determined by LAL gel clotting assay.
Isotype	IgG2ak
Storage Instruction	Store at 2-8°C. Do not freeze this product.

TARGET INFORMATION

Gene ID	17067 100041546
Gene Symbol	Ly6c1 Ly6c2
Uniprot ID	LY6C1_MOUSE LY6C2_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	