

Anti-CD45.2/Ptprc antibody [104.2] (STJA0001762)

STJA0001762

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Receptor-type tyrosine-protein phosphatase C is suitable for use in Flow Cytometry and Immunohistochemistry research applications.
Applications	FC/IHC
Host/Source	Mouse
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	104.2
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300
Isotype	IgG2ak
Storage	Store at 2-8°C. Do not freeze this product.
Instruction	

TARGET INFORMATION

Gene ID	19264
Gene Symbol	Ptprc
Uniprot ID	PTPRC_MOUSE
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	