

**Anti-TCR Beta antibody [H57-597 (HB218)] (STJA0001768)**

STJA0001768

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Armenian Hamster monoclonal antibody anti-TCR Beta is suitable for use in Universal, Flow Cytometry, Immunohistochemistry and Immunoprecipitation research applications.
<b>Applications</b>	Depletion/FC/IHC/IP
<b>Host/Source</b>	Armenian Hamster
<b>Reactivity</b>	Mouse

**PRODUCT PROPERTIES**

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	H57-597 (HB218)
<b>Concentration</b>	0.5 mg/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	0.2 µm filtered in PBS, pH7.2. Azide Free (AF) /Low Endotoxin (LE) : Contains no stabilizers or preservatives. Endotoxin level is < 2 EU/mg as determined by LAL gel clotting assay.
<b>Isotype</b>	IgG
<b>Storage</b>	Store at 2-8°C. Do not freeze this product.
<b>Instruction</b>	

**TARGET INFORMATION**

<b>Gene ID</b>
<b>Gene Symbol</b>
<b>Uniprot ID</b>
<b>Immunogen</b>
<b>Immunogen</b>
<b>Region</b>
<b>Specificity</b>
<b>Immunogen</b>
<b>Sequence</b>