

**Anti-CD154/CD40LG antibody [5C8] {Biotin} (STJA0001808)**

STJA0001808

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-CD40 ligand is suitable for use in Flow Cytometry research applications.
<b>Applications</b>	FC
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human/Monkey

**PRODUCT PROPERTIES**

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	5C8
<b>Concentration</b>	0.5 mg/mL
<b>Conjugation</b>	Biotin
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
<b>Isotype</b>	IgG2ak
<b>Storage Instruction</b>	Store at 2-8°C. Do not freeze this product.

**TARGET INFORMATION**

<b>Gene ID</b>	<a href="#">959</a>
<b>Gene Symbol</b>	<a href="#">CD40LG</a>
<b>Uniprot ID</b>	<a href="#">CD40L_HUMAN</a>
<b>Immunogen</b>	
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	