

Anti-HLA-A, B, C antibody [W6/32] {Biotin} (STJA0001813)
STJA0001813**GENERAL INFORMATION**

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-HLA-A, B, C is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	W6/32
Concentration	0.5 mg/mL
Conjugation	Biotin
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2ak
Storage Instruction	Store at 2-8°C. Do not freeze this product.

TARGET INFORMATION

Gene ID	3106
Gene Symbol	HLA-B
Uniprot ID	HLAB_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	