

**Anti-Cd274/PD-L1 antibody [10F.9G2] (STJA0001820)**

STJA0001820

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rat monoclonal antibody anti-Programmed cell death 1 ligand 1 is suitable for use in Flow Cytometry, Immunohistochemistry, Immunofluorescence and Universal research applications.
<b>Applications</b>	FC/IHC/IF/Block
<b>Host/Source</b>	Rat
<b>Reactivity</b>	Mouse

**PRODUCT PROPERTIES**

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	10F.9G2
<b>Concentration</b>	0.5 mg/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	0.2 µm filtered in PBS, pH7.2. Azide Free (AF) /Low Endotoxin (LE) : Contains no stabilizers or preservatives. Endotoxin level is < 2 EU/mg as determined by LAL gel clotting assay.
<b>Isotype</b>	IgG2bk
<b>Storage</b>	Store at 2-8°C. Do not freeze this product.
<b>Instruction</b>	

**TARGET INFORMATION**

<b>Gene ID</b>	60533
<b>Gene Symbol</b>	Cd274
<b>Uniprot ID</b>	PD1L1_MOUSE
<b>Immunogen</b>	
<b>Immunogen</b>	
<b>Region</b>	
<b>Specificity</b>	
<b>Immunogen</b>	
<b>Sequence</b>	