

Anti-Cd274/PD-L1 antibody [10F.9G2] (STJA0001821)

STJA0001821

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Primary antibodies |
| Short Description | Rat monoclonal antibody anti-Programmed cell death 1 ligand 1 is suitable for use in Flow Cytometry, Immunohistochemistry and Immunofluorescence research applications. |
| Applications | FC/IHC/IF |
| Host/Source | Rat |
| Reactivity | Mouse |

PRODUCT PROPERTIES

| | |
|-----------------------|---|
| Clonality | Monoclonal |
| Clone ID | 10F.9G2 |
| Concentration | 0.5 mg/mL |
| Conjugation | Unconjugated |
| Purification | |
| Dilution Range | |
| Formulation | PBS with 0.05% Proclin300 |
| Isotype | IgG2bk |
| Storage | Store at 2-8°C. Do not freeze this product. |
| Instruction | |

TARGET INFORMATION

| | |
|--------------------|-------------|
| Gene ID | 60533 |
| Gene Symbol | Cd274 |
| Uniprot ID | PD1L1_MOUSE |
| Immunogen | |
| Immunogen | |
| Region | |
| Specificity | |
| Immunogen | |
| Sequence | |