

Anti-CD49d/ITGA4 antibody [9F10] {FITC} (STJA0001861)

STJA0001861

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Integrin alpha-4 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	9F10
Concentration	5 Mu L/Test
Conjugation	FITC
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG1k
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	3676
Gene Symbol	ITGA4
Uniprot ID	ITA4_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	