

Anti-CD11b/ITGAM antibody [ICRF44] {PerCP/Cyanine5.5} (STJA0001878)

STJA0001878

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Integrin alpha-M is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	ICRF44
Concentration	5 Mu L/Test
Conjugation	PerCP/Cyanine5.5
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG1k
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	3684
Gene Symbol	ITGAM
Uniprot ID	ITAM_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	