

Anti-CD134/TNFRSF4 antibody [Ber-ACT35] {FITC} (STJA0001895)
STJA0001895**GENERAL INFORMATION**

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Tumor necrosis factor receptor superfamily member 4 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	Ber-ACT35
Concentration	5 Mu L/Test
Conjugation	FITC
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG1k
Storage	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.
Instruction	

TARGET INFORMATION

Gene ID	7293
Gene Symbol	TNFRSF4
Uniprot ID	TNR4_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	