

Anti-CD303/CLEC4C antibody [201A] {APC} (STJA0001926)

STJA0001926

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-C-type lectin domain family 4 member C is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	201A
Concentration	5 Mu L/Test
Conjugation	APC
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2ak
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	170482
Gene Symbol	CLEC4C
Uniprot ID	CLC4C_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen	
Sequence	