

Anti-CD39/ENTPD1 antibody [A1] {PE/Cyanine5} (STJA0001931)

STJA0001931

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Ectonucleoside triphosphate diphosphohydrolase 1 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	A1
Concentration	5 Mu L/Test
Conjugation	PE/Cyanine5
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG1k
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	953
Gene Symbol	ENTPD1
Uniprot ID	ENTP1_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	