

Anti-CD90/THY1 antibody [5E10] {PE/Cyanine5.5} (STJA0001940)

STJA0001940

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Thy-1 membrane glycoprotein is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	5E10
Concentration	5 µg L/Test
Conjugation	PE/Cyanine5.5
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG1k
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	7070
Gene Symbol	THY1
Uniprot ID	THY1_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	