

Anti-IL17A antibody [BL168] {PE} (STJA0001960)
STJA0001960**GENERAL INFORMATION**

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Interleukin-17A is suitable for use in IConfocal MicroscopyC research applications.
Applications	ICFC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	BL168
Concentration	5 Mu L/Test
Conjugation	PE
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG1k
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	3605
Gene Symbol	IL17A
Uniprot ID	IL17_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	