

Anti-CD273/PDCD1LG2 antibody [24F.10C12] {APC} (STJA0001965)

STJA0001965

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Programmed cell death 1 ligand 2 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	24F.10C12
Concentration	5 Mu L/Test
Conjugation	APC
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2ak
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	80380
Gene Symbol	PDCD1LG2
Uniprot ID	PD1L2_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen	
Sequence	