

Anti-Cd14 antibody [Sa14-2] {AF647} (STJA0001975)

STJA0001975

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Monocyte differentiation antigen CD14 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	Sa14-2
Concentration	0.5 mg/mL
Conjugation	AF647
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2ak
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	12475
Gene Symbol	Cd14
Uniprot ID	CD14_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen	
Sequence	