

Anti-CD150/Slamf1 antibody [TC15-12F12.2] {FITC} (STJA0001976)

STJA0001976

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Signaling lymphocytic activation molecule is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	TC15-12F12.2
Concentration	5 Mu L/Test
Conjugation	FITC
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2a, lambda
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	27218
Gene Symbol	Slamf1
Uniprot ID	SLAF1_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen	
Sequence	