

Anti-CD335/Ncr1 antibody [29A1.4] {PE/Cyanine7} (STJA0002009)
STJA0002009**GENERAL INFORMATION**

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Natural cytotoxicity triggering receptor 1 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	29A1.4
Concentration	5 Mu L/Test
Conjugation	PE/Cyanine7
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2ak
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	17086
Gene Symbol	Ncr1
Uniprot ID	NCTR1_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen	
Sequence	