

Anti-Fcεr1α antibody [MAR-1] {APC} (STJA0002044)

STJA0002044

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Armenian Hamster monoclonal antibody anti-High affinity immunoglobulin epsilon receptor subunit alpha is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Armenian Hamster
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	MAR-1
Concentration	5 Mu L/Test
Conjugation	APC
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG
Storage	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.
Instruction	

TARGET INFORMATION

Gene ID	14125
Gene Symbol	Fcεr1α
Uniprot ID	FCERA_MOUSE
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	