

Anti-Ly6a antibody [D7] {AF488} (STJA0002092)

STJA0002092

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Lymphocyte antigen 6A-2/6E-1 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	D7
Concentration	0.5 mg/mL
Conjugation	AF488
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2ak
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	110454
Gene Symbol	Ly6a
Uniprot ID	LY6A_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	