

Anti-CD366/HAVCR2 antibody [RMT3-23] {APC} (STJA0002097)

STJA0002097

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Hepatitis A virus cellular receptor 2 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	RMT3-23
Concentration	5 Mu L/Test
Conjugation	APC
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2ak
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	84868
Gene Symbol	HAVCR2
Uniprot ID	HAVR2_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen	
Sequence	