

Anti-Cd38 antibody [NIMR5] (STJA0002101)

STJA0002101

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-ADP-ribosyl cyclase/cyclic ADP-ribose hydrolase 1 is suitable for use in Flow Cytometry and Immunohistochemistry research applications.
Applications	FC/IHC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	NIMR5
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300
Isotype	IgG2ak
Storage	Store at 2-8°C. Do not freeze this product.
Instruction	

TARGET INFORMATION

Gene ID	12494
Gene Symbol	Cd38
Uniprot ID	CD38_MOUSE
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	