

**Anti-IL4 antibody [11B11] (STJA0002157)**

STJA0002157

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rat monoclonal antibody anti-Interleukin-4 is suitable for use in Universal, Universal, IConfocal MicroscopyC, Western Blot, Immunohistochemistry and Immunoprecipitation research applications.
<b>Applications</b>	Neut/Stim/ICFC/WB/IHC/IP
<b>Host/Source</b>	Rat
<b>Reactivity</b>	Mouse

**PRODUCT PROPERTIES**

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	11B11
<b>Concentration</b>	0.5 mg/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	0.2 µm filtered in PBS, pH7.2. Azide Free (AF) /Low Endotoxin (LE) : Contains no stabilizers or preservatives. Endotoxin level is < 2 EU/mg as determined by LAL gel clotting assay.
<b>Isotype</b>	IgG1k
<b>Storage</b>	Store at 2-8°C. Do not freeze this product.
<b>Instruction</b>	

**TARGET INFORMATION**

<b>Gene ID</b>	<a href="#">16189</a>
<b>Gene Symbol</b>	<a href="#">IL4</a>
<b>Uniprot ID</b>	<a href="#">IL4_MOUSE</a>
<b>Immunogen</b>	
<b>Immunogen</b>	
<b>Region</b>	
<b>Specificity</b>	
<b>Immunogen</b>	
<b>Sequence</b>	