

Anti-IL-5 antibody [TRFK5] (STJA0002163)

STJA0002163

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-IL-5 is suitable for use in Universal, Universal, IConfocal MicroscopyC, Western Blot and Immunohistochemistry research applications.
Applications	Depletion/Neut/ICFC/WB/IHC
Host/Source	Rat
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	TRFK5
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	0.2 µm filtered in PBS, pH7.2. Azide Free (AF) /Low Endotoxin (LE) : Contains no stabilizers or preservatives. Endotoxin level is < 2 EU/mg as determined by LAL gel clotting assay.
Isotype	IgG1k
Storage	Store at 2-8°C. Do not freeze this product.
Instruction	

TARGET INFORMATION

Gene ID	3567 16191
Gene Symbol	IL5 IL5
Uniprot ID	IL5_HUMAN IL5_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	