

Anti-IL-5 antibody [TRFK5] (STJA0002164)

STJA0002164

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-IL-5 is suitable for use in IConfocal MicroscopyC research applications.
Applications	ICFC
Host/Source	Rat
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	TRFK5
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300
Isotype	IgG1k
Storage Instruction	Store at 2-8°C. Do not freeze this product.

TARGET INFORMATION

Gene ID	16191 3567
Gene Symbol	IL5 IL5
Uniprot ID	IL5_MOUSE IL5_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	