

Anti-IL6 antibody [MQ2-13A5] {PE} (STJA0002170)

STJA0002170

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Primary antibodies |
| Short Description | Rat monoclonal antibody anti-Interleukin-6 is suitable for use in IConfocal MicroscopyC research applications. |
| Applications | ICFC |
| Host/Source | Rat |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|----------------------------|--|
| Clonality | Monoclonal |
| Clone ID | MQ2-13A5 |
| Concentration | 5 Mu L/Test |
| Conjugation | PE |
| Purification | |
| Dilution Range | |
| Formulation | Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant. |
| Isotype | IgG1k |
| Storage Instruction | Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product. |

TARGET INFORMATION

| | |
|---------------------------|---------------------------|
| Gene ID | 3569 |
| Gene Symbol | IL6 |
| Uniprot ID | IL6_HUMAN |
| Immunogen | |
| Immunogen Region | |
| Specificity | |
| Immunogen Sequence | |