

**Anti-CD44 antibody [Hermes-1] (STJA0002208)**

STJA0002208

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rat monoclonal antibody anti-CD44 antigen is suitable for use in Flow Cytometry and Immunofluorescence research applications.
<b>Applications</b>	FC/IF
<b>Host/Source</b>	Rat
<b>Reactivity</b>	Human

**PRODUCT PROPERTIES**

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	Hermes-1
<b>Concentration</b>	0.5 mg/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	PBS with 0.05% Proclin300
<b>Isotype</b>	IgG2ak
<b>Storage Instruction</b>	Store at 2-8°C. Do not freeze this product.

**TARGET INFORMATION**

<b>Gene ID</b>	960
<b>Gene Symbol</b>	CD44
<b>Uniprot ID</b>	CD44_HUMAN
<b>Immunogen</b>	
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen</b>	
<b>Sequence</b>	