

Anti-CD44 antibody [Hermes-1] {AF647} (STJA0002213)

STJA0002213

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-CD44 antigen is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	Hermes-1
Concentration	5 Mu L/Test
Conjugation	AF647
Purification	
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2ak
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	960
Gene Symbol	CD44
Uniprot ID	CD44_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	