

Anti-H-2 antibody [M1/42] {PerCP/Cyanine5.5} (STJA0002226)
STJA0002226**GENERAL INFORMATION**

| | |
|--------------------------|-----------------------------------------------------------------------------------------------|
| Product Type | Primary antibodies |
| Short Description | Rat monoclonal antibody anti-H-2 is suitable for use in Flow Cytometry research applications. |
| Applications | FC |
| Host/Source | Rat |
| Reactivity | Mouse |

PRODUCT PROPERTIES

| | |
|----------------------------|--------------------------------------------------------------------------------------------------------------|
| Clonality | Monoclonal |
| Clone ID | M1/42 |
| Concentration | 0.2 mg/mL |
| Conjugation | PerCP/Cyanine5.5 |
| Purification | |
| Dilution Range | |
| Formulation | Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant. |
| Isotype | IgG2ak |
| Storage Instruction | Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product. |

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen Region
Specificity
Immunogen Sequence