

Anti-CD3E antibody [UCHT1] (STJA0002232)

STJA0002232

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-T-cell surface glycoprotein CD3 epsilon chain is suitable for use in Universal, Universal, Universal, Flow Cytometry, Immunohistochemistry and Western Blot research applications.
Applications	Activ/Depletion/Inhibition/FC/IHC/WB
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	UCHT1
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	0.2 µm filtered in PBS, pH7.2. Azide Free (AF) /Low Endotoxin (LE) : Contains no stabilizers or preservatives. Endotoxin level is < 2 EU/mg (determined by LAL gel clotting assay).
Isotype	IgG1k
Storage	Store at 2-8°C. Do not freeze this product.
Instruction	

TARGET INFORMATION

Gene ID	916
Gene Symbol	CD3E
Uniprot ID	CD3E_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	