

Anti-Foxp3 antibody [3G3] {FITC} (STJA0002249)

STJA0002249

GENERAL INFORMATION


| | |
|--------------------------|--|
| Product Type | Primary antibodies |
| Short Description | Mouse monoclonal antibody anti-FoXp3 is suitable for use in IConfocal MicroscopyC research applications. |
| Applications | ICFC |
| Host/Source | Mouse |
| Reactivity | Mouse |

PRODUCT PROPERTIES

| | |
|----------------------------|--|
| Clonality | Monoclonal |
| Clone ID | 3G3 |
| Concentration | 5 Mu L/Test |
| Conjugation | FITC |
| Purification | |
| Dilution Range | |
| Formulation | Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant. |
| Isotype | IgG1k |
| Storage Instruction | Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product. |

TARGET INFORMATION

| | |
|---------------------------|-----------------------------|
| Gene ID | 50943 |
| | 20371 |
| Gene Symbol | FOXP3 |
| | Foxp3 |
| Uniprot ID | FOXP3_HUMAN |
| | FOXP3_MOUSE |
| Immunogen | |
| Immunogen Region | |
| Specificity | |
| Immunogen Sequence | |



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081